

PARAFFIN WAX WORKS

QUICK HEAT Therapeutic Paraffin Kit



PARAFFIN WAX WORKS

Item #301476
tng worldwide • New Hudson, MI 48165
paraffinwaxworks.com • 800.520.6000
©2014 TNG • Made in China

IMPORTANT INSTRUCTIONS — RETAIN FOR FUTURE USE

Congratulations on your purchase of the **Paraffin Wax Works Therapeutic Quick Heat Paraffin Kit (Unit)**. Please follow the simple instructions and precautions before operating the unit. The unit comes with a one-year limited warranty. If for any reason you have questions or issues, please contact us toll-free 800-520-6000.

Instructions For Use:

1. Please place the Unit on safe and flat surface before using it.
2. Plug the Unit into a grounded outlet.
3. Add 3 pounds of Paraffin Wax Works paraffin wax to Unit. The Unit will hold a total of 5 pounds.
4. Turn the switch to Hi. Allow up to 45 minutes for full melting of wax to take place. Once wax is melted you can change to Low if you have heat sensitivity.
5. Clean hands thoroughly. For first time users, gently dip one finger into paraffin to experience the warmth sensation. Then gently dip hand into unit keeping fingers straight, palm first. Remove hand immediately and allow excess wax to drip off. Repeat twice.
6. Immediately wrap hand with Paraffin Wax Works plastic liner. Then cover plastic liner with thermal mitts. Repeat with other hand. Note: Can also be used on feet. Can also be used on elbows. When using on body, always use a brush to apply paraffin to skin.
7. Allow up to 20 minutes for paraffin to moisturize the skin. This will help blood circulation, increase oxygen supply to the blood and help stimulate the vitality of cells improving metabolism.
8. After 20 minutes, remove thermal mitt. Then remove plastic liner starting at the wrist and gently pull down off the hand taking as much paraffin wax as possible. For best results do not reuse paraffin wax.
9. After treatment, turn off the power and pull plug from outlet. Wax will harden and be good for next use.
10. For optimal results, do treatments at least weekly.

Precautions:

1. Do not use outdoors or near water.
2. Keep out of reach of children and pets.
3. Always keep paraffin wax at minimum level to prevent burning.
4. For best results, refill Unit with Paraffin Wax Works paraffin beads.
5. Do not melt paraffin wax in a microwave oven.

Specifications:

Model #	301476
Volume	5 pounds
Wattage	180W
Voltage	110V/60Hz
Temperature	158° F

Limited Warranty

Paraffin Wax Works offers a one-year limited warranty for defects only. Not responsible for Unit due to negligence, dropping, outdoor use or operating near water. If you experience any issues with the Unit please contact us at 800-520-6000. Do not contact or return to the retailer that you purchased from. In case of warranty replacement, please have your proof of purchase available.